

Product Information	
Catalog No.	EL10212
Product name	Anti-Atezolizumab Neutralizing Antibody ELISA kit
Applications	ELISA
Accession	Q9NZQ7
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	15.625-1, 000 ng/ml
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	MPDL3280A, RG7446, 1380723-44-3
Shipping	2-8 deg C
Note	For Research Use Only.