

Product Information	
Catalog No.	EL10250
Product name	Anti-Nivolumab ELISA Kit
Applications	ELISA
Accession	Q15116
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Alternative Names	BMS-936558, MDX-1106, ONO-4538, 946414-94-4
Shipping	2-8 deg C
Note	For Research Use Only.