

Product Information	
Catalog No.	EL10373
Product name	Insulin Lispro ELISA Kit
Applications	ELISA
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	33.12 ng/mL
Range	78.125 - 5, 000 ng/mL
Recovery	80-120%
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Shipping	2-8 deg C
Note	For Research Use Only.