

Product Information	
Catalog No.	PA13740
Product name	Anti-Mouse CD117/KIT Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Mast/stem cell growth factor receptor Kit (SCFR) (EC 2.7.10.1) (Proto-oncogene c-Kit) (Tyrosine-protein kinase Kit) (CD antigen CD117)
Immunogen	E. coli - derived recombinant Mouse CD117/KIT (Ser38-Gly513).
Accession	P05532
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.