


Product Information

Catalog No.	AT0044
Product name	Anti-APP4-tag Antibody (SAA1990)
Applications	ELISA, SPR
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1990
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	APP, Amyloid beta A4 protein, Alzheimer disease amyloid protein, Beta-amyloid precursor protein, Protease nexin-II
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-APP4-tag Antibody (SAA1990)

SDS-PAGE