

Product Information	
Catalog No.	AT0066
Product name	Anti-EGFP Monoclonal Antibody (1A190), HRP
Applications	ELISA, IF, WB
Species reactivity	Water jellyfish
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A190
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Purification	Affinity chromatography.
Immunogen	Recombinant EGFP Protein.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Alternative Names	EGFP
Note	For research use only.