


Product Information	
Catalog No.	BS10048
Product name	Anti-Human beta-NGF Reference Antibody (Fulranumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Beta-nerve growth factor, NGFB, Beta-NGF, NGF
Expression system	Mammalian Cells
Accession	P01138
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Alternative Names	4D4, AMG-403, JNJ-42160443, 902141-80-4
Species	Human
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Fulranumab.

SDS-PAGE