


Product Information	
Catalog No.	BS10145
Product name	Anti-Human P-Selectin/CD62P Reference Antibody (Inclacumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.68mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Leukocyte-endothelial cell adhesion molecule 3, Platelet activation dependent granule-external membrane protein, Granule membrane protein 140, GMP-140, CD62P, SELP, GRMP, PADGEM, CD62 antigen-like family member P, LECAM3, P-selectin, GMRP
Expression system	Mammalian Cells
Accession	P16109
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	LC1004-002, RO4905417, 1256258-86-2
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Torapsel (42-89-glycoprotein; WAY164339) is a fusion protein with immunoglobulin. Torapsel can be used to research the prevention of ischemia reperfusion injury.

Image Data



SDS-PAGE for Research Grade Inclacumab.

SDS-PAGE