


Product Information


Catalog No.	BS10151
Product name	Anti-Human CEACAM-5/CD66e Reference Antibody (Cergutuzumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-Kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Meconium antigen 100, CEACAM5, CEA, CD66e, Carcinoembryonic antigen-related cell adhesion molecule 5, Carcinoembryonic antigen
Expression system	Mammalian Cells
Accession	P06731
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	CEA-IL-2v, RG7813, RO-6895882, 1509916-03-3
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data

Cergutuzumab binds with CD66e/CEA/CEACAM5

Detects Human CD66e/CEA/CEACAM5 in indirect ELISAs.





SDS-PAGE for Research Grade Cergutuzumab.

SDS-PAGE