


Product Information	
Catalog No.	BS10316
Product name	Anti-Human IgE/IGHE Reference Antibody (Talizumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-nd
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ig epsilon chain C region, IGHE, Ig epsilon chain C region ND, Immunoglobulin heavy constant epsilon
Expression system	Mammalian Cells
Accession	P01854
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	C21/AL-901, HU-901, Hu901, TNX-901, 380610-22-0
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SEC-HPLC detection for Research Grade Quilizumab.

Bioactivity