


Product Information	
Catalog No.	BS10321
Product name	Anti-Human Complement C5 Reference Antibody (Crovalimab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CPAMD4, C5, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 4, Complement C5
Expression system	Mammalian Cells
Accession	P01031
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	RG-6107, RO7112689, SKY-59, 1917321-26-6
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Omoprubart is an anti-C5 human?IgG2-G4 κ monoclonal antibody.

Image Data



SEC-HPLC detection for Research Grade Crovalimab.

Bioactivity