


Product Information	
Catalog No.	BS10400
Product name	Anti-Human FOLH1/PSMA & CD3E Bispecific Reference Antibody (Pasotuxizumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	(scFv-heavy-kappa)-(scFv-heavy-lambda)
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	NAALAD1, Prostate-specific membrane antigen, NAALADase I, PSMA, PSM, Pteroylpoly-gamma-glutamate carboxypeptidase, Glutamate carboxypeptidase 2, N-acetylated-alpha-linked acidic dipeptidase I, Glutamate carboxypeptidase II, FOLH1, FGCP, FOLH, Folylpoly-gamma-glutamate carboxypeptidase, Folate hydrolase 1, GCPII, Membrane glutamate carboxypeptidase, Cell growth-inhibiting gene 27 protein, mGCP, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain
Expression system	Mammalian Cells
Accession	Q04609 & P07766
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	Bispecific, BAY2010112, BITEMT112, MT-112, 1442657-12-6
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Voxalatamab (JNJ-63898081) is a bispecific IgG4 antibody that binds PSMA and CD3 and has anti-cancer activity for use in prostate cancer research.

Image Data

Pasotuxizumab binds with CD3E



Bioactivity

Detects CD3E in indirect ELISAs.