


Product Information


Catalog No.	BS10513
Product name	Anti-HRSV F/F protein Reference Antibody (HNK20)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	HRSV-A
Host species	Mouse
Clonality	Monoclonal
Isotype	IgA
Clone ID	HNK20
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Fusion glycoprotein F2, p27, Intervening segment, Pep27, Peptide 27, Fusion glycoprotein F1
Expression system	Mammalian Cells
Accession	P03420
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



Bioactivity

SEC-HPLC detection for Research Grade Anti-HRSV F/F protein Antibody (HNK20).



SDS-PAGE for Research Grade Anti-HRSV F/F protein Antibody (HNK20).

SDS-PAGE