


Product Information	
Catalog No.	BS10542
Product name	Anti-Human CXCR5 Reference Antibody (Sar113244, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG4-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CXC-R5, CD185, MDR-15, C-X-C chemokine receptor type 5, CXCR5, CXCR-5, MDR15, Monocyte-derived receptor 15, Burkitt lymphoma receptor 1, BLR1
Expression system	Mammalian Cells
Accession	P32302
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	SAR113244
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade SAR113244

SDS-PAGE