



Product Information	
Catalog No.	BS10633
Product name	Anti-Human KLKB1/CLK3 Reference Antibody (Lanadelumab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	3.16mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Plasma prekallikrein, KLK3, Fletcher factor, KLKB1, Kininogenin, Plasma kallikrein, PKK
Expression system	Mammalian Cells
Accession	P03952
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	X124-G01, DX-2930, 1426055-14-2
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



Bioactivity

SEC-HPLC detection for Research Grade Lanadelumab.



SDS-PAGE for Research Grade Lanadelumab.

SDS-PAGE