


Product Information	
Catalog No.	BS10659
Product name	Anti-Human CD19 Reference Antibody (Loncastuximab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Chimeric
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.94 mM Citricacid, 16.7mM Trisodiumcitratetrihydrate, 6.2% sucrose, pH 6.5, 0.01% Tween80
Concentration	4.55mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD19, Differentiation antigen CD19, T-cell surface antigen Leu-12, B-lymphocyte antigen CD19, B-lymphocyte surface antigen B4
Expression system	Mammalian Cells
Accession	P15391
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	ADCT-402, RB4v1.2-SG-3249, 1875032-68-0
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SEC-HPLC detection for Research Grade Loncastuximab.

Bioactivity