

Product Information	
Catalog No.	BS10663
Product name	Anti-Human HLA-G Reference Antibody (HLA-G-0104)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HLA-G, HLA-G, sHLA-G, HLA G antigen, MHC class I antigen G, HLA-6.0, HLA class I histocompatibility antigen, alpha chain G
Expression system	Mammalian Cells
Accession	P17693
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	HLA-G-0104, HLA-G-0031
Note	For research use only. Not suitable for clinical or therapeutic use.