


Product Information	
Catalog No.	BS10680
Product name	Anti-Human ErbB3/Her3 Reference Antibody (Barecetamab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-Lambda
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	ERBB3, HER3, Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor HER3, Receptor tyrosine-protein kinase erbB-3
Expression system	Mammalian Cells
Accession	P21860
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	ISU-104, ISU104, 2275727-74-5
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Nezutatug is an anti-HER3 antibody, and used for cancer research.

Image Data



SDS-PAGE for Research Grade Eluvixtamab

SDS-PAGE