


Product Information	
Catalog No.	BS10700
Product name	Anti-Human Proprotein Convertase 9/PCSK9 Reference Antibody (Recaticimab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-Kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Proprotein convertase subtilisin/kexin type 9, Subtilisin/kexin-like protease PC9, NARC-1, PC9, Proprotein convertase 9, Neural apoptosis-regulated convertase 1, NARC1, PCSK9
Expression system	Mammalian Cells
Accession	Q8NBP7
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	SHR-1209, 2361290-85-7
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Frovocimab (LY3015014) (INN) is a monoclonal antibody designed for the treatment of hypercholesterolemia. This drug is being developed by Eli Lilly and Company.

Image Data



SDS-PAGE for Research Grade Recaticimab.

SDS-PAGE