



Product Information	
Catalog No.	BS10702
Product name	Anti-Human IL-33 Reference Antibody (Tozorakimab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-Lambda
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	NF-HEV, IL-1F11, C9orf26, IL-33, Interleukin-1 family member 11, IL33, Nuclear factor from high endothelial venules, NFHEV, Interleukin-33, IL1F11
Expression system	Mammalian Cells
Accession	O95760
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	MEDI-3506, 2376858-66-9
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Tozorakimab.

SDS-PAGE



SEC-HPLC detection for Research Grade Tozorakimab.

Bioactivity