

Product Information	
Catalog No.	BS10703
Product name	Anti-Human CD33 & CD3E Bispecific Reference Antibody (Eluvixtamab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	(scFv-heavy-kappa)-(scFv-heavy-lambda)
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Myeloid cell surface antigen CD33, Sialic acid-binding Ig-like lectin 3, Siglec-3, CD33, gp67, SIGLEC3, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain
Expression system	Mammalian Cells
Accession	P20138 & P07766
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	Bispecific, AMG-330, MT-114, 1679391-73-1
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Eluvixtamab is an anti-CD3E/CD33 monoclonal antibody.

Image Data



SDS-PAGE for Research Grade Eluvixtamab

SDS-PAGE