

Product Information	
Catalog No.	BS10738
Product name	Anti-Human APOC3 Reference Antibody (Stt-5058, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Humanized
Clonality	Monoclonal
Isotype	IgG1-Kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	APOC3, Apolipoprotein C3, Apo-CIII, Apolipoprotein C-III, ApoC-III
Expression system	Mammalian Cells
Accession	P02656
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	ARGX-116
Note	For research use only. Not suitable for clinical or therapeutic use.