


Product Information	
Catalog No.	BS10912
Product name	Anti-Human MICB&MICA Reference Antibody (Iph4301, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MICA, MIC-A, MHC class I polypeptide-related sequence A, PERB11.1, MHC class I polypeptide-related sequence B, MIC-B, MICB, PERB11.2
Expression system	Mammalian Cells
Accession	Q29983, Q29980
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	IPH4301
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade IPH4301.

SDS-PAGE