


Product Information	
Catalog No.	BS10971
Product name	Anti-Human CD276/B7-H3 Reference Antibody (DS-5573a)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	DS-5573a
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	4.3mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD276, B7-H3, Costimulatory molecule, B7H3, 4Ig-B7-H3, B7 homolog 3, CD276 antigen
Expression system	Mammalian Cells
Accession	Q5ZPR3
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	DS-5573a, DS5573a, DS5573a
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CD276/B7-H3 (DS-5573a).

SDS-PAGE