


| Product Information   |                                                                                                                                                                                                                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | BS11096                                                                                                                                                                                                                                                                                           |
| Product name          | Anti-Human CD38/ADPRC 1 Reference Antibody (SAR442085)                                                                                                                                                                                                                                            |
| Applications          | ELISA, Bioactivity: FACS, Functional assay, Research in vivo                                                                                                                                                                                                                                      |
| Species reactivity    | Human                                                                                                                                                                                                                                                                                             |
| Host species          | Human                                                                                                                                                                                                                                                                                             |
| Clonality             | Monoclonal                                                                                                                                                                                                                                                                                        |
| Isotype               | IgG                                                                                                                                                                                                                                                                                               |
| Clone ID              | SAR442085                                                                                                                                                                                                                                                                                         |
| Form                  | Liquid                                                                                                                                                                                                                                                                                            |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                                                                                                                                                                |
| Concentration         | 1mg/ml                                                                                                                                                                                                                                                                                            |
| Purity                | >95% purity as determined by SDS-PAGE.                                                                                                                                                                                                                                                            |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                                                                                                                                                               |
| Target                | Cyclic ADP-ribose hydrolase 1, ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1, cADPr hydrolase 1, 2'-phospho-ADP-ribosyl cyclase/2'-phospho-cyclic-ADP-ribose transferase, T10, 2'-phospho-cyclic-ADP-ribose transferase, CD38, ADPRC 1, ADP-ribosyl cyclase 1, 2'-phospho-ADP-ribosyl cyclase |
| Expression system     | Mammalian Cells                                                                                                                                                                                                                                                                                   |
| Accession             | P28907                                                                                                                                                                                                                                                                                            |
| Endotoxin level       | Please contact the lab for this information.                                                                                                                                                                                                                                                      |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.                                                                                                   |
| Alternative Names     | SAR442085                                                                                                                                                                                                                                                                                         |
| Note                  | For research use only. Not suitable for clinical or therapeutic use.                                                                                                                                                                                                                              |

### Image Data



SDS-PAGE for Research Grade Anti-Human CD38/ADPRC 1 (SAR442085)

SDS-PAGE