

Product Information

Catalog No.	BS11214
Product name	Research Grade Acazicolcept

Background

Host species:

Isotype: Fusion-[ICOSLGN-terminalfragment(1-122)]2-IGHG1Fc(Fragmentconstant)

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality:

Applications: Research Grade Biosimilar

Target: CD28,TP44,T-cell-specificsurfaceglycoproteinCD28,CD278,ICOS,InducibleT-cellcostimulator,Activation-induciblelymphocyteimmunomediatorymolecule,AILIM

Purification: Protein A/G purified from cell culture supernatant.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Alternative Names: ALPN-101,ICOSLvIgD-Fc,CAS:2270247-50-0

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



Detects CD28 in indirect ELISAs.



SDS-PAGE for Research Grade Acazicolcept

SDS-PAGE