

Product Information

Catalog No.	BS11216
Product name	Research Grade Asunercept

Background

Host species:

Isotype: Fusion-[FAS(Fascellsurfacedeathreceptor,TNFRSF6,tumornecrosisfactorreceptor(TNFR)superfamilymember6,FAS1,APO-1,CD95)]2-IGHG1Fc(Fragmentconstant)

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 1mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality:

Applications: Research Grade Biosimilar

Target: Receptor-bindingFasLectodomain,APL,APT1LG1,FasLICD,Fasantigenligand,CD178,sFasL,CD95L,APTL,TNFSF6,Apoptosisantigenligand,SPA,FASLG,FasL,SPPL2A-processedFasLform,Tumornecrosisfactorligandsuperfamilymember6,SolubleFasligand,Fasligand,CD95ligand,CD95-L,FASL

Purification: Protein A/G purified from cell culture supernatant.


Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: APG-101,CAS:1450882-18-4

Note: For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CX43/GJA1 (ALMB-0168)

SDS-PAGE