


Product Information	
Catalog No.	BS11261
Product name	Anti-Human CXCL13/BCA-1/BLC Reference Antibody (TJX7)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	TJX7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1
Expression system	Mammalian Cells
Accession	O43927
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	TJX7
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CXCL13/BCA-1/BLC (TJX7)

SDS-PAGE