

Product Information	
Catalog No.	BS11268
Product name	Anti-Human CD184/CXCR4 Reference Antibody (PF-06747143)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	PF-06747143
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.11mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Stromal cell-derived factor 1 receptor, LCR1, HM89, LPS-associated protein 3, LAP-3, CXC-R4, CXCR-4, SDF-1 receptor, CD184, NPYRL, LESTR, Fusin, FB22, C-X-C chemokine receptor type 4, CXCR4, Lipopolysaccharide-associated protein 3, Leukocyte-derived seven transmembrane domain receptor
Expression system	Mammalian Cells
Accession	P61073
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	PF-06747143, PF06747143, PF06747143
Note	For research use only. Not suitable for clinical or therapeutic use.