

Product Information	
Catalog No.	BS11269
Product name	Anti-Human CD185/CXCR5 Reference Antibody (BAY-924)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	BAY-924
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.24mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CXC-R5, CD185, MDR-15, C-X-C chemokine receptor type 5, CXCR5, CXCR-5, MDR15, Monocyte-derived receptor 15, Burkitt lymphoma receptor 1, BLR1
Expression system	Mammalian Cells
Accession	P32302
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	BAY-924, BAY924, BAY924
Note	For research use only. Not suitable for clinical or therapeutic use.