


Product Information	
Catalog No.	BS11317
Product name	Anti-Human F11 Reference Antibody (BAY 1831865)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	BAY 1831865
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Plasma thromboplastin antecedent, PTA, Coagulation factor XI, F11, FXI
Expression system	Mammalian Cells
Accession	P03951
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	BAY-1831865, BAY1831865, BAY1831865, FXIa, FactorXI
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human F11 (BAY 1831865)

SDS-PAGE