

Product Information	
Catalog No.	BS11345
Product name	Anti-Human LGALS9 Reference Antibody (LYT-200)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4-kappa
Clone ID	LYT-200
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	4.53mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ecalectin, LGALS9, Galectin-9, Gal-9, Tumor antigen HOM-HD-21
Expression system	Mammalian Cells
Accession	O00182
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	LYT-200, LYT200, LYT200
Note	For research use only. Not suitable for clinical or therapeutic use.