


Product Information	
Catalog No.	BS11352
Product name	Anti-RSV F/Fusion glycoprotein F0 Reference Antibody (R4.C6)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human respiratory syncytial virus A (strain A2)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	R4.C6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F0, Fusion glycoprotein F2, p27, Intervening segment, Pep27, Peptide 27, Fusion glycoprotein F1
Expression system	Mammalian Cells
Accession	P03420
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data




SDS-PAGE for Research Grade Anti-RSV F/Fusion glycoprotein F0 (R4.C6)

SDS-PAGE

Anti-RSV F/Fusion glycoprotein F0 (R4.C6)
binds with Fusion glycoprotein F0

Detects recombinant HRSV-A2 Pre-F/Fusion glycoprotein F0 (Catalog No: VK421011) in indirect ELISA



Bioactivity