


Product Information	
Catalog No.	BS11376
Product name	Anti-Human CD20 Reference Antibody (Blontuvetmab, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Dog
Host species	Caninized
Clonality	Monoclonal
Isotype	IgG2-kappa-lambda
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD20/MS4A1
Expression system	Mammalian Cells
Accession	Q3C2E2
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	AT-004, VT-007, 1608112-78-2
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Blontuvetmab.

SDS-PAGE