


Product Information	
Catalog No.	BS11410
Product name	Anti-Human TGF-beta 1 Reference Antibody (Linavonkibart, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transforming growth factor beta-1 proprotein, LAP, TGF-beta-1, TGFB1, TGFB
Expression system	Mammalian Cells
Accession	P01137
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2640981-30-0
Note	For research use only. Not suitable for clinical or therapeutic use.

Background

Linavonkibart (SRK 181) is a high-affinity, fully human antibody that selectively binds to latent TGFβ1 and inhibits its activation. Linavonkibart plays an important role in cancer.

Image Data



SDS-PAGE for Research Grade Linavonkibart

SDS-PAGE