


Product Information	
Catalog No.	BS11417
Product name	Anti-Human PD-1&TIGIT Reference Antibody (Rilvegostomig, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa/lambda
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	VSTM3, TIGIT, V-set and transmembrane domain-containing protein 3, VSIG9, V-set and immunoglobulin domain-containing protein 9, T-cell immunoreceptor with Ig and ITIM domains, Programmed cell death protein 1, Protein PD-1, hPD-1, PD1, PDCD1, CD279
Expression system	Mammalian Cells
Accession	Q495A1 & Q15116
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	2640305-01-5
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Rilvegostomig

SDS-PAGE