

Product Information

Catalog No.	BS11494
Product name	Research Grade Balugrastim

Background

Host species: Human

Isotype:

Fusion-ALB(albumin,human serum albumin,HSA)25-609-CSF3(colony stimulating factor 3,G-CSF,granulocyte colony-stimulating factor,G-CSF)

Species reactivity: Human

Immunogen:

Form: Liquid

Storage instructions:

Storage buffer: 0.01M PBS, pH7.4.

Concentration: 0.5mg/ml

Purity: >95% as determined by SDS-PAGE.

Clonality:

Applications: Research Grade Biosimilar

Target: G-CSF receptor, CSF3R, GCSFR, G-CSF-R, CD114, Granulocyte colony-stimulating factor receptor

Purification: Purified by Ion Exchange Chromatography.

Endotoxin level: Expression system: Mammalian Cells

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Alternative Names: ALB-rG-CSF, G-CSF-HAS, HSA-rG-CSF, EGRANLI[®], Φ , NEUGRANIN[®], Φ , ALBUGRANIN[®], Φ

Note: For research use only. Not suitable for clinical or therapeutic use.