


Product Information	
Catalog No.	BS11503
Product name	Anti-Human B7-H4 Reference Antibody (Xmt-1660, RUO)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	V-set domain-containing T-cell activation inhibitor 1, B7h.5, B7 homolog 4, B7-H4, T-cell costimulatory molecule B7x, B7H4, VTCN1, Protein B7S1, Immune costimulatory protein B7-H4
Expression system	Mammalian Cells
Accession	Q7Z7D3
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade XMT-1660.

SDS-PAGE