



Product Information	
Catalog No.	BS11525
Product name	Anti-Human GFRAL Reference Antibody (NGM120)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	NGM120
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	C6orf144, GDNF family receptor alpha-like, GFRAL
Expression system	Mammalian Cells
Accession	Q6UXV0
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	NGM120
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human GFRAL (NGM120).

SDS-PAGE



Western blot analysis was performed using anti-GFRAL monoclonal antibody at 1ug/ml on various samples.

Lane 1: recombinant human GFRAL (Catalog No: HS720011), 50ng

Lane 2: recombinant human GFRAL (Catalog No: HS720011), 30ng

Western blot