

Product Information	
Catalog No.	BS11530
Product name	Anti-Human CD357/TNFRSF18/GITR Reference Antibody (JNJ-64164711)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	JNJ-64164711
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Activation-inducible TNFR family receptor, GITR, TNFRSF18, CD357, Tumor necrosis factor receptor superfamily member 18, AITR, Glucocorticoid-induced TNFR-related protein
Expression system	Mammalian Cells
Accession	Q9Y5U5
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	JNJ-64164711, JNJ64164711, JNJ64164711
Note	For research use only. Not suitable for clinical or therapeutic use.