


Product Information	
Catalog No.	BS11619
Product name	Anti-Human IFNAR1 Reference Antibody (QX 006N)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	QX 006N
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Interferon alpha/beta receptor 1, Cytokine receptor family 2 member 1, Cytokine receptor class-II member 1, CRF2-1, IFN-R-1, IFNAR1, IFNAR, IFN-alpha/beta receptor 1, Type I interferon receptor 1
Expression system	Mammalian Cells
Accession	P17181
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	QX-006N, QX006N, QX006N
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human IFNAR1 (QX 006N)

SDS-PAGE