


Product Information	
Catalog No.	BS11656
Product name	Anti-Human IL1RAP Reference Antibody (CSC012)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	CSC012
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-1R3, C3orf13, Interleukin-1 receptor 3, Interleukin-1 receptor accessory protein, IL1RAP, IL-1R-3, IL-1 receptor accessory protein, IL1R3, IL-1RAcP
Expression system	Mammalian Cells
Accession	Q9NPH3
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	CSC012
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human IL1RAP (CSC012)

SDS-PAGE