


Product Information	
Catalog No.	BS11667
Product name	Anti-Human CD360/IL21R Reference Antibody (ATR-107)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human, Mouse
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-lambda
Clone ID	ATR-107
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL21R, NILR, Interleukin-21 receptor, IL-21R, IL-21 receptor, CD360, Novel interleukin receptor
Expression system	Mammalian Cells
Accession	Q9HBE5
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	ATR-107, ATR107, ATR107
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data

Detects Recombinant Human CD360/IL21R Protein, C-His (Cat No.: HV076011) in indirect ELISA.



Detects Recombinant Mouse CD360/IL21R Protein, C-His (Cat No.: MV076011) in indirect ELISA.

