



Product Information	
Catalog No.	BS11730
Product name	Anti-Human CD305/LAIR1 Reference Antibody (IO 106)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1-kappa
Clone ID	IO 106
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.83mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Leukocyte-associated immunoglobulin-like receptor 1, LAIR-1, hLAIR1, LAIR1, CD305
Expression system	Mammalian Cells
Accession	Q6GTX8
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	IO-106, IO106, IO106
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CD305/LAIR1 (IO 106).

SDS-PAGE




Western blot analysis was performed using anti-CD305/LAIR1 monoclonal antibody at 1ug/mL on various samples.

Lane 1: recombinant human CD305/LAIR1 (Catalog No: HV765011)

Lane 2: negative control

Western blot

**Anti-Human CD305/LAIR1 (IO 106)
binds with CD305**



Detects recombinant CD305 (Catalog No: HV765011) in indirect ELISA.