


Product Information	
Catalog No.	BS11747
Product name	Anti-Human CSF1/M-CSF Reference Antibody (PD-0360324)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	PD-0360324
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CSF1, Macrophage colony-stimulating factor 1, CSF-1, Lanimostim, M-CSF, MCSF
Expression system	Mammalian Cells
Accession	P09603
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	PD-0360324, PD0360324, PD0360324
Note	For research use only. Not suitable for clinical or therapeutic use.


Image Data



SDS-PAGE for Research Grade Anti-Human CSF1/M-CSF (PD-0360324)

SDS-PAGE

Anti-Human CSF1/M-CSF (PD-0360324)
binds with CSF1/M-CSF



Detects recombinant human CSF1/M-CSF (Catalog No: HY515011) in indirect ELISA