


Product Information	
Catalog No.	BS11772
Product name	Anti-Human NONO Reference Antibody (PAT-LM1)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	PAT-LM1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	54 kDa nuclear RNA- and DNA-binding protein, NRB54, p54nrb, NMT55, NonO protein, p54(nrb), 55 kDa nuclear protein, NONO, DNA-binding p52/p100 complex, 52 kDa subunit, Non-POU domain-containing octamer-binding protein
Expression system	Mammalian Cells
Accession	Q15233
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	PAT-LM1, PATLM1, PATLM1
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human NONO (PAT-LM1)

SDS-PAGE