


Product Information	
Catalog No.	BS11875
Product name	Anti-Human SIGLEC15/CD33L3 Reference Antibody (BSI-060T)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	BSI-060T
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	5.11mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD33L3, Siglec-15, Sialic acid-binding Ig-like lectin 15, SIGLEC15, CD33 antigen-like 3
Expression system	Mammalian Cells
Accession	Q6ZMC9
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	BSI-060T, BSI060T, BSI060T
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human SIGLEC15/CD33L3 (BSI-060T).

SDS-PAGE