


Product Information	
Catalog No.	BS11918
Product name	Anti-Human TGFB1/TGF-beta-1 Reference Antibody (SAR-439459)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	SAR-439459
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transforming growth factor beta-1 proprotein, LAP, TGF-beta-1, TGFB1, TGFB
Expression system	Mammalian Cells
Accession	P01137
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	SAR-439459, SAR439459, SAR439459
Note	For research use only. Not suitable for clinical or therapeutic use.


Image Data



SDS-PAGE for Research Grade Anti-Human TGFB1/TGF-beta-1 (SAR-439459)

SDS-PAGE

Anti-Human TGFB1/TGF-beta-1 (SAR-439459) binds with TGFB1



Detects recombinant TGFB1 (Catalog No: HF977011) in indirect ELISA.