



Product Information	
Catalog No.	BS11920
Product name	Anti-Human TGFB1/TGF-beta-1 Reference Antibody (SRK-181)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4
Clone ID	SRK-181
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.35mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transforming growth factor beta-1 proprotein, LAP, TGF-beta-1, TGFB1, TGFB
Expression system	Mammalian Cells
Accession	P01137
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	SRK-181, SRK181, SRK181
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data

**Anti-Human TGFB1/TGF-beta-1 (SRK-181)
binds with TGFB1**

Detects Recombinant Human TGFB1/TGF-beta-1 Protein, C-His (Cat No.: HF977011) in indirect ELISA.



SDS-PAGE for Research Grade Anti-Human TGFB1/TGF-beta-1 (SRK-181).

SDS-PAGE