


Product Information	
Catalog No.	BS11952
Product name	Anti-Human CD120b/TNFRSF1B/TNFR2 Reference Antibody (BC011)
Applications	ELISA, Bioactivity: FACS, Functional assay, Research in vivo
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG
Clone ID	BC011
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1mg/ml
Purity	>95% purity as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TNF-R2, TNFR-II, p80 TNF-alpha receptor, Etanercept, TNFBR, TBP-2, TBPII, Tumor necrosis factor receptor superfamily member 1B, CD120b, p75, TNFRSF1B, Tumor necrosis factor receptor 2, TNF-RII, TNFR2, Tumor necrosis factor receptor type II
Expression system	Mammalian Cells
Accession	P20333
Endotoxin level	Please contact the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C for short-term storage (1-2 weeks). Store at -20°C for up to 12 months. For long-term storage, store at -80°C.
Alternative Names	BC011
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for Research Grade Anti-Human CD120b/TNFRSF1B/TNFR2 (BC011)

SDS-PAGE